


# Computed Tomography and Stationary Digital Tomosynthesis Pulmonary Vessel Alignment

3D to 2D to 3D registration

Natalie Johnston  
Summer Internship 2022



# Background

Benefits of Tomosynthesis over CT include greater practical accessibility and less radiation exposure for the patient


CT to tomosynthesis vessel registration can enhance guided medical procedures by increasing availability and decreasing the patient's exposure to radiation

Driving clinical problem: image-guided lung biopsy


Registering pre-operative 3D CT with intra-operative 3D Tomosynthesis reconstruction

# Data


- CT volume
- Tomosynthesis projection images
- Tomosynthesis reconstruction volume
- Geometries of 29 emitter positions




CT Volume [3D]



Segmented Vessels [3D]



Tomosynthesis  
Projection\* [2D]



Tomosynthesis  
Reconstruction [3D]

\*29 images produced by 29 emitter geometries

# Tomosynthesis Projections





Image: Shan et.al 2015

# Goal: register CT 3D vessels with 2D tomosynthesis projection image




# Method

1. Get CT source points from vessel extraction
2. Get tomosynthesis reconstruction vessel points from hand-drawn annotation
3. Make mask with drawn tomosynthesis reconstruction points
4. Apply the mask to each of the tomosynthesis projection images
5. Get masked tomosynthesis images (from 29 emitter positions)
6. Project and register CT source points with the mask applied to the tomosynthesis projection



# 1. CT vessel centerline points from extraction

- itk-TubeTK
  - Vessel segmentation
  - Tube spatial objects
- Spatial object position in object space becomes source points for vessel registration



## 2. TomoRecon points from hand drawn annotation

- ImageViewer
  - Annotate prominent vessels in the tomosynthesis reconstruction 2D slices
- Points extracted from the overlay are used as a mask



3. Make mask with drawn tomoRecon points
4. Apply the mask to each of the tomosynthesis projection images
5. Get masked tomosynthesis images (from 29 emitter positions)

- Using the known emitter positions (geometries), use the points from the overlay to create 29 masks for the 29 tomosynthesis projections



\*tomo projection images: 02 and 28

## 6. Project and register CT source points with the mask applied to the tomosynthesis projection (algorithm)

Small perturbation of values in x,  
repeat until x tolerance is satisfied


```
leastSquaresOptimizer:  
    for each point:  
        for each emitterPosition:  
            transform point given x  
            project point onto tomoProjMask  
                using corresponding emitter position  
            store voxel value at  
                coordinates of  
                projected point  
  
            avg voxel value of point across emitter positions  
            avg voxel value of points  
  
maximize voxel values (optimize to brightest parts)
```

```
return x=[xTranslation, yTranslation, zTranslation, zRotation, yRotation, xRotation]
```


**Inputs:**  
x=[0, 0, 0, 0, 0, 0]  
u=sourcePoints  
y=emitterGeometry  
mask\_list  
x\_scale  
size

## 6. Project and register CT source points with the mask applied to the tomosynthesis projection (algorithm)


LeastSquaresOptimizer



# Results




# Verification: Combined Image




Segmented vessels from CT


Transform  
with solution  
from  
optimization




Transformed vessels from CT and tomosynthesis  
reconstruction




# Verification: Combined Image Slices




Pre-Registration



Post-Registration




Pre-Registration




Post-Registration

Patient #2

# Verification: Combined Image Slices




Pre-Registration




Post-Registration

Patient #1

# Verification: Combined Image Slices




Pre-Registration




Post-Registration

Patient #3

# Verification: Combined Image Slices





Pre-Registration



Post-Registration

Patient #4

# Verification: Combined Image Slices



Pre-Registration


Post-Registration

Patient #5

# Acknowledgements



Dr. Stephen Aylward



Dr. Yueh Lee